



Analyse Exomes CLOVES


Dr. Thomas Steimlé


Mars 2023

Laboratoire d'oncohématologie de l'hôpital Necker

Materials & Methods

Alias	Cases	Patho	Normal
AUZAL	--	612200191032	612200128560
BASPA	--	2007No80798	2007No80327
BOUDE-1	--	2006No80622	2006No80641
BOUDE-2	--	2006No80656	2006No80641
BRIAN	--	2105N190483	2105N170599
DEVCL	--	2010N260761	2010N260754
DUMEL	--	612200588717	612200588556
FARMO	--	2105N210788	2105N210782
GOSJU	--	612120459843	612120454105
KRAMA	--	1911No40786	1911No40784
LABPH	--	2006No80633	2006No80640
LASPH	--	612120335709	612120316267
LECFA	--	2008N170574	2007N270441
LHYET	--	2006No80768	2006No80767
MASFA	--	612100070595	612100070580
MOKRA	--	1811N280725	1501N300498
MULAY	--	6612112860	61210092424
OUATI	--	2104N160342	2104N130615
PLANY	--	2103N110343	2103N110124
REMMA	--	2104N160304	2104N150297
ROQRO-1	--	2006No80629	2006No80644
ROQRO-2	--	612200572392	6122005723638
SAMAR	--	1712No40796	612120165753
VERPY	--	612120515645	612120504794

 Séquençage selon la méthode *Agilent SureSelect Human All Exon V7 panel* sur automates *Illumina Next/NovaSeq*.

 Analyse bioinformatique:

Alignement (hg19 bwa mem v0.7.17-r1188) et déduplication UMI (umi_tools v1.1.1)

Callers :



mutect2 (GATK v4.2)



strelka (Illumina v2.9.10)




lancet (NY Genome Center v1.1.0)

En suivant les modes opératoires disponibles (cf. diapos supplémentaires).

Results


↔ Alignment – proportions of targeted sequences at given depth by case




Case		< 100x (%)			
—	LASPH_N	100	—	MOKRA_P	7
—	LASPH_P	100	—	OUATI_N	7
—	BRIAN_P	41	—	REMMAN	7
—	BOUDE-1_P	26	—	REMMAN_P	6
—	KRAMA_N	26	—	LABPH_P	6
—	KRAMA_P	20	—	PLANY_N	6
—	MULAY_P	19	—	ROQRO-1_N	6
—	MULAY_N	18	—	MOKRA_N	6
—	DEVCL_N	15	—	ROQRO-1_P	5
—	MASFA_N	11	—	BASPA_N	5
—	DEVCL_P	11	—	GOSJU_N	5
—	BOUDE-2_P	11	—	ROQRO-2_P	5
—	LECFA_N	11	—	VERPY_P	5
—	BASPA_P	11	—	BOUDE-1_N	5
—	SAMAR_P	10	—	BOUDE-2_N	5
—	LECFA_P	10	—	SAMAR_N	5
—	BRIAN_N	10	—	VERPY_N	5
—	LHYET_N	10	—	DUMEL_P	5
—	FARMO_N	9	—	DUMEL_N	4
—	LHYET_P	9	—	ROQRO-2_N	4
—	MASFA_P	8	—	GOSJU_P	4
—	OUATI_P	7	—	AUZAL_N	4
—	PLANY_P	7	—	LABPH_N	3
—	FARMO_P	7	—	AUZAL_P	3




Variant calling



Cases on the left with an abnormal mutiounal burden and LASPH which is has been sequenced as a panel have been filtered out.




Filtering – Proportion of filtered calls for each case






Filtering – Rules and resulting allelic fraction distribution

Number of calls filtered by each five applied filters



Number of calls for each allelic fraction after filtering






Results – Recurrent genes known as tumor suppressors or proto-oncogenes.

Alias	Gene	Variant	VAF (%)
BASPA	ALK	p.Gly1202GlufsTer56	2.9
REMMA	ALK	p.Lys1267Asn	1.4
REMMA	ARID1A	p.Thr290Pro	9.3
PLANY	ARID1A	p.Arg1223His	1
MULAY	ARID1A	p.Met1759Ile	2.4
MULAY	ARID1A	p.Ala532Val	5.9
PLANY	KMT2C	p.His4278Asn	1.4
FARMO	KMT2C	p.Ala117Thr	1
DEVCL	KMT2C	intron	0.6
BASPA	NCOA3	p.Gln758Ter	2.5
MOKRA	NCOA3	p.Pro115Leu	1.1
REMMA	NCOR1	p.Arg627Leu	1.7
MULAY	NCOR1	intron	3.7
MULAY	NCOR1	p.Gly1340Cys	1.3
GOSJU	NOTCH1	p.Gln1837Ter	2.7
GOSJU	NOTCH1	p.Gly251Asp	1.2
ROQRO-2	NOTCH1	p.Ala1418Thr	1.1
BOUDE-2	NOTCH1	intron	1.4
ROQRO-2	NOTCH2	p.Gln2006Ter	1.2
BASPA	NOTCH2	p.His1208Asn	1.5
MOKRA	PIK3CA	p.Glu545Lys	3.5
DEVCL	PIK3CA	p.Glu542Lys	4.2
SAMAR	SMARCA4	p.Met797Ile	2
MASFA	SMARCA4	p.Phe939Leu	1.9






Results – Over-representation analysis – KEGG pathways | WikiPathways | Reactome








Results – Over-representation analysis – KEGG – Focal Adhesion






Complexes:

2.7.1.149
1-phosphatidylinositol-5-phosphate 4-kinase
2.7.1.153
phosphatidylinositol-4,5-bisphosphate 3-kinase catalytic subunit alpha/beta/delta
2.7.1.154
phosphatidylinositol-4-phosphate 3-kinase
2.7.1.67
phosphatidylinositol 4-kinase A
2.7.1.68
1-phosphatidylinositol-4-phosphate 5-kinase
3.1.3.-
phosphatidylinositol 4-phosphatase
3.1.3.25
myo-inositol-1(or 4)-monophosphatase
3.1.3.36
inositol polyphosphate 5-phosphatase INPP5B/F
3.1.3.48/3.1.3.64/3.1.3.95
myotubularin-related protein 3/4
3.1.3.56
inositol polyphosphate 5-phosphatase INPP5J/K
3.1.3.64/3.1.3.95
myotubularin-related protein 1/2
3.1.3.66
inositol polyphosphate-4-phosphatase
3.1.3.86
phosphatidylinositol-3,4,5-trisphosphate 5-phosphatase 1
3.1.4.11
phosphatidylinositol phospholipase C, beta





Results – Over-representation analysis – KEGG – Choline metabolism in cancer






Results – Confirmation panel




Confirmation panel results: 3646 positions have been included in the confirmation panel, 65 have been confirmed.



Results – Confirmation panel



The exome has a better depth for the variant positions we wanted to confirm



The sequencing depths at positions we wanted to confirm are deeper in the exome. In the confirmation panel, confirmed calls positions have significantly deeper depth.

Conclusions

- ⇒ 17 exomes analysés, environ 300 mutations par cas.
- ⇒ Mutations retrouvées majoritairement à des VAF allant de 1 à 2 %
- ⇒ Récurrence de mutations sur les gènes connus pour être impliqués dans le cancer: *ALK*, *ARID1A*, *KMT2C*, *NCOA3*, *NCOR1*, *NOTCH1*, *NOTCH2*, *PIK3CA* et *SMARCA4*.
- ⇒ Net enrichissement de mutations touchant des gènes impliqués dans le métabolisme des Inositol Phosphates.
- ⇒ Panel de confirmation insuffisamment couvert cependant il confirme 65 mutations.

Supplementaries



Source code : <https://github.com/nygenome/lancet>

```
1 lancet --ref hg19.fa
2     --tumor {TUMOR.BAM}
3     --normal {NORMAL.BAM}
4     --bed {AGILENT_REGIONS_V7.BED}
5     --num-threads 31
6     > {OUTPUT.VCF}
```

Listing 1: lancet – bash version

Manta: <https://github.com/Illumina/manta>

```
1 configManta.py --exome
2   --referenceFasta hg19.fa
3   --tumorBam {TUMOR.BAM}
4   --normalBam {NORMAL.BAM}
5   --callRegions {AGILENT_REGIONS_V7.BED}
6   --runDir /tmp/...
```

Listing 2: manta – bash version

Strelka: <https://github.com/Illumina/strelka>

```
1 configureStrelkaSomaticWorkflow.py --exome
2   --referenceFasta hg19.fa
3   --tumorBam {TUMOR.BAM}
4   --normalBam {NORMAL.BAM}
5   --callRegions {AGILENT_REGIONS_V7.BED}
6   --indelCandidates {MANTA_candidateSmallIndels.vcf.gz}
7   --runDir /tmp/...
```

Listing 3: strelka – bash version



From : <https://gatk.broadinstitute.org/hc/en-us/articles/360035531132>

```
1 gatk-4.2.1.0/gatk --java-options -Xmx32g
2   Mutect2
3   -R hg19.fa
4   -L {AGILENT_REGIONS_V7.BED}
5   -I {TUMOR.BAM}
6   -I {NORMAL.BAM}
7   --normal-sample {NORMAL}
8   --germline-resource af-only-gnomad.raw.sites.chr.vcf
9   --panel-of-normals pon.vcf.gz
10  --f1r2-tar-gz {TUMOR_tumoral_f1r2.tar.gz}
11  -O {OUTPUT.VCF}
```

Listing 4: Mutect2 – seule étape adaptée – bash version